

PATIENT LIFTING SLING





FAST

Our large warehouse and sophisticated fulfilment systems mean we can deliver with speed. This is supported by our online customer services, making it easy to track orders and re-order.

INTELLIGENT

Years of research and development have allowed us to find care solutions that are affordable, innovative and useful. Our products cover manual handling pressure care and patient safety.

CARE

We know our products inside and out and we want to share this knowledge ans expertise. We help our customers find the right cre solutions, as well as offering ongoing training, support and product maintenance.

TABLE OF CONTENTS


Loop Slings Access Hygiene Loop Quickfit GP Loop Quickfit Deluxe Leg Loop Full Back GP Head Support	4
Loop Slings Loop Universal MultiFit	8
Clip Slings Clip Comfort Head Support Clip Comfort Access Head Support	10
Active Slings Loop Stand Deluxe Loop Stand Deluxe Transport Deluxe Stand-Aid Comfort Clip	12
FIM Score Assessment Tool	14
Notes	15

LOOP SLINGS

Oxford® Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

LIFT TYPES	SPREADER BAR	SLING SWL
Advance Presence Midi 180	Loop	227 kg



Loop Sling Access Hygiene - Padded	Size	Enable Code	Α	В	С	D
	XS	LG6245	600mm	260mm	570mm	N/A
8	S	LG6246	700mm	280mm	660mm	N/A
	М	LG6247	850mm	310mm	890mm	N/A
	L	LG6248	1050mm	330mm	1000mm	N/A
	XL	LG6249	1300mm	390mm	1160mm	N/A

Loop Access Sling Hygiene Head Support - Padded	Size	Enable Code	Α	В	С	D
	XS	LG6235	590mm	270mm	1180mm	440mm
	S	LG6236	690mm	300mm	1360mm	490mm
	М	LG6237	830mm	310mm	1550mm	540mm
	L	LG6238	1060mm	280mm	1600mm	590mm
	XL	LG6239	1260mm	340mm	1860mm	670mm

Loop Sling Quickfit General Purpose - Padded	Size	Enable Code	A	В	С	D
	XS	Special	640mm	470mm	960mm	N/A
<u> </u>	S	LG6047	730mm	560mm	1120mm	N/A
6	М	LG6059	900mm	620mm	1210mm	N/A
	L	LG6071	1110mm	700mm	1400mm	N/A
1	XL	LG6079	1390mm	910mm	1780mm	N/A

Loop Sling Quickfit General Purpose - Mesh	Size	Enable Code	А	В	С	D
	XS	LG6538	640mm	470mm	960mm	N/A
- A B	S	LG6051	730mm	560mm	1120mm	N/A
	М	LG6063	900mm	620mm	1210mm	N/A
	L	LG6075	1110mm	700mm	1400mm	N/A
1	XL	LG6540	1390mm	910mm	1780mm	N/A

Loop Sling Quickfit General Purpose - Polyester	Size	Enable Code	Α	В	С	D
	XS	Special	640mm	470mm	960mm	N/A
B	S	LG6043	730mm	560mm	1120mm	N/A
	М	LG6055	900mm	620mm	1210mm	N/A
	L	LG6067	1110mm	700mm	1400mm	N/A
	XL	Special	1390mm	910mm	1780mm	N/A

SLING	SIZE	COLOR		ı	RECO	MMEN	IDED	WEIG	HT RA	ANGE	(kg)		FIM
SLING	SIZE	COLOR	15	34	45	57	68	80	90	125	136	227	SCORE
Loop Quickfit Deluxe Leg	S	Red			LG	5010							
Material: Polyester	М	Yellow					LG!	5014					1 & 2
	L	Green							LG!	5018			









*All slings on this page may be used for amputees

Loop Sling Quickfit Deluxe Leg - Polyester	Size	Enable Code	Α	В	С	D
	S	LG5010	770mm	630mm	920mm	N/A
	М	LG5014	890mm	770mm	1120mm	N/A
	L	LG5018	1030mm	850mm	1180mm	N/A

Loop Sling Quickfit Deluxe Leg - Mesh	Size	Enable Code	Α	В	С	D
	S	LG5012	770mm	630mm	920mm	N/A
3	М	LG5016	890mm	770mm	1120mm	N/A
	L	LG5020	1030mm	850mm	1180mm	N/A

Loop Sing Quickfit Deluxe Leg - Padded	Size	Enable Code	Α	В	С	D
	S	LG6095	770mm	630mm	920mm	N/A
	М	LG6101	890mm	770mm	1120mm	N/A
	L	LG6107	1030mm	850mm	1180mm	N/A

Loop Sling Full Back GP Head Support - Polyester	Size	Enable Code	А	В	С	D
1	XS	LG6117	810mm	860mm	1310mm	N/A
	S	LG5030	980mm	890mm	1370mm	N/A
	М	LG5032	1180mm	1060mm	1580mm	N/A
	L	LG5034	1360mm	1160mm	1760mm	N/A
\ /	XL	LG6125	1690mm	1340mm	2090mm	N/A

Loop Sling Full Back GP Head Support - Mesh	Size	Enable Code	А	В	С	D
\/	XS	LG6127	810mm	860mm	1310mm	N/A
	S	LG6129	980mm	890mm	1370mm	N/A
	М	LG6131	1180mm	1060mm	1580mm	N/A
Į V	L	LG6133	1360mm	1160mm	1760mm	N/A
\ /	XL	LG6135	1690mm	1340mm	2090mm	N/A

LOOP SLINGS


Oxford® Universal Loop Sling MultiFit Colour/Size Chart - Standard Size

Universal sling suitable for toilet, general transfers and amputees. Recommended weight range should only be used as a guideline as other factors can affect sizing.

LIFT TYPES	SPREADER BAR	SLING SWL
Presence Midi 180 Advance	Loop	227 kg

Sizing Guide: Oxford® MultiFit Slings

		1	2	3	4	5	6	7	8	9	10	11	12
ion	A	540mm	600mm	660mm	720mm	780mm	840mm	910mm	980mm	1040mm	1110mm	1200mm	1270mm
nensi	В	480mm	520mm	570mm	620mm	670mm	750mm	810mm	840mm	880mm	910mm	930mm	950mm
Din	С	800mm	880mm	970mm	1070mm	1170mm	1310mm	1410mm	1460mm	1530mm	1570mm	1600mm	1630mm















Sizing Guide: Oxford® Universal Loop Sling MultiFit

	Height (cm)	Girth (cm)	Max SWL	Size	Colour Code
	75 - 95	45 - 55	227 kg / 35st	1 (XS)	LG6263
Paediatric	95 - 110	50 - 60	227 kg / 35st	2 (XS)	LG6264
(XS - S)	110 - 130	55 - 70	227 kg / 35st	3 (S)	LG6265
	130 - 150	65 - 80	227 kg / 35st	4 (S)	LG6266
	135 - 160	70 - 85	227 kg / 35st	5 (M)	LG6267
	145 - 170	75 - 90	227 kg / 35st	6 (M)	LG6268
	155 - 180	85 - 100	227 kg / 35st	7 (M)	LG6269
Adult	160 - 190	95 - 115	227 kg / 35st	8 (M)	LG6270
Male (M - XL)	160 - 190	105 - 130	227 kg / 35st	9 (L)	LG6271
	160 - 190	125 - 150	227 kg / 35st	10 (L)	LG6272
	160 - 190	145 - 175	284 kg / 44st	11 (XL)	LG6273
	160 - 190	170 - 200	284 kg / 44st	12 (XL)	LG6274
	130 - 155	85 - 100	227 kg / 35st	5 (M)	LG6267
	135 - 160	95 - 110	227 kg / 35st	6 (M)	LG6268
Adult	140 - 165	105 - 125	227 kg / 35st	7 (M)	LG6269
Female	140 - 165	120 - 140	227 kg / 35st	8 (M)	LG6270
(M - XL)	140 - 165	135 - 155	227 kg / 35st	9 (L)	LG6271
	140 - 165	150 - 175	227 kg / 35st	10 (L)	LG6272
	140 - 165	170 - 200	284 kg / 44st	11 (XL)	LG6273

Measuring the Client

You will need to take two measurements to find the most suitable sling size for your client; height and girth. See guidance below whilst referring to the diagram and table above for assistance.

1. Height

Take your client's full height as shown. It may prove more accurate to measure height from a recumbent (lying) position.

2. Girth

For female clients, measure around the hips. For male clients and children, measure around the chest just under the arms.

Universal MultiFit

This Universal MultiFit incorporates over-sized leg supports to specifically cater for clients with larger than average hips/thighs. Should you find a client's girth size is too large to match the height in the table above, try these for comparison.

CLIP SLINGS

Oxford® Clip Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

LIFT TYPES	CRADLE	SLING SWL
Presence Advance (Shallow Cradle)	Clip	227 kg



Clip Sling Comfort Head Support - Polyester	Size	Enable Code	Α	В	С	D
	XS	LG6372	810mm	860mm	1310mm	N/A
	S	LG6373	980mm	890mm	1370mm	N/A
	М	LG6374	1180mm	1060mm	1580mm	N/A
	L	LG6375	1360mm	1160mm	1760mm	N/A
.	XL	LG6376	1690mm	1340mm	2090mm	N/A

Clip Sling Comfort Head Support - Mesh	Size	Enable Code	Α	В	С	D
	XS	LG6382	810mm	860mm	1310mm	N/A
B	S	LG6383	980mm	890mm	1370mm	N/A
	М	LG6384	1180mm	1060mm	1580mm	N/A
	L	LG6385	1360mm	1160mm	1760mm	N/A
Į. Ų Ų	XL	LG6386	1690mm	1340mm	2090mm	N/A

Clip Sling Comfort Head Support - Padded	Size	Enable Code	А	В	С	D
1	XS	LG6377	810mm	860mm	1310mm	N/A
<u>A</u>	S	LG6378	980mm	890mm	1370mm	N/A
	М	LG6379	1180mm	1060mm	1580mm	N/A
	L	LG6380	1360mm	1160mm	1760mm	N/A
Į V	XL	LG6381	1690mm	1340mm	2090mm	N/A

Clip Sling Comfort Access Head Support - Mesh	Size	Enable Code	A	В	С	D
	S	LG6394	690mm	300mm	1360mm	490mm
	М	LG6395	830mm	310mm	1550mm	540mm
Ų Ų.	L	LG6396	1060mm	280mm	1600mm	590mm

ACTIVE SLINGS

Oxford® Loop and Clip Style Sling Color/Size Chart - Standard Size

Recommended weight range should only be used as a guideline as other factors can affect sizing.

LIFT TYPES	CRADLE	SLING SWL
Ascend Journey	Loop 😘	200 kg

SLING	CLZE	COLOR		F	RECO	1MEN	IDED \	WEIGI	HT RA	NGE (kg)		FIM
SLING	SIZE	COLOR	15	34	45	57	68	80	90	125	136	200	SCORE
Loop Stand Deluxe	S	Red			LS	5112							
Material: Polyester	М	Yellow					LS	5114				-	704
8))	L	Green							LS	5110			3 & 4
	XL	Blue		•	***************************************	••••		•	*************		LS5110	6	
Loop Stand Deluxe Transport	S	Red		L	G618	5		•••••	•••••	•	***************************************	•••••	
Material: Polyester	М	Yellow					LG6	5184					2, 3 & 4
	L	Green							LG	6186			

LIFT TYPES	CRADLE	SLING SWL
Oxford Switch Oxford Up	Clip	200 kg

SLING	SIZE	COLOR	RECOMMENDED WEIGHT RANGE (kg)								FIM		
			15	34	45	57	68	80	90	125	136	200	SCORE
Deluxe Stand-Aid Comfort Clip	S	Red			LS	5118							
Material: Polyester	М	Yellow					LS5	120					7.0.4
	L	Green	LS5122					3 & 4					
	XL	Blue			***************************************		•			l	LS512	4	

Loop Sling Stand Deluxe - Polyester	Size	Enable Code	Α	В	С	D
	S	LS5112	720mm	310mm	740mm	N/A
	М	LS5114	830mm	330mm	900mm	N/A
	L	LS5110	1070mm	340mm	1220mm	N/A
	XL	LS5116	1250mm	410mm	1400mm	N/A

Loop Sling Stand Deluxe Transport - Polyester	Size	Enable Code	A	В	С	D
	S	LG6185	670mm	290mm	860mm	N/A
	М	LG6184	840mm	340mm	1010mm	N/A
	L	LG6186	1100mm	400mm	1060mm	N/A

Deluxe Stand-Up / Standing / Transport Sling, Loop, Clip	Size	Enable Code	Α	В	С	D
	S	LS5118	720mm	310mm	740mm	N/A
	М	LS5120	830mm	330mm	900mm	N/A
	L	LS5122	1070mm	340mm	1220mm	N/A
	XL	LS5124	1250mm	410mm	1400mm	N/A

FIM SCORE EXPLAINED

FIM (Functional Independence Measure) is a widely used tool, particularly within the field of rehabilitation and involves the classification of a person's mobility according to 7 levels of function, from independence to total assistance. (Granger CV, Hamilton BB, Linacare JM, Heinemann AW, Wright BD, Performance profiles of the functional independence measure. Am J Phys Med Rehabill 1993; 72:84-9).

Independent

Another person is not required for the activity. (NO HELPER)

7. Complete Independence

All of the tasks described as making up the activity are typically performed safely, without modification, assistive devices or aids and within a reasonable time.

6. Modified Independence

Activity requires one or more than one of the following: an assistive device, more than reasonable time, or safety (risk) considerations

Dependent

Another person is required for either supervision or physical assistance in order for the activity to be performed (REQUIRES HELPER).

5. Supervision or set-up

Person requires no more help than standby, cueing or coaxing without physical contact, or helper sets up needed items or applies orthoses.

4. Minimal Assistance

With physical contact, the person requires no more help than touching and person expends 75% or more of the effort.

3. Moderate Assistance

Person requires more help than touching alone or expends half (50%) or more (up to 75%) of the effort.

Complete Dependence

The person expends less than half (50%) of the effort. Maximal or total assistance is required, or the activity is not performed.

2. Maximal Assistance

Subject expends less than 50% of the effort, but at least 25%.

1. Total Assistance

Subject expends less than 25% of the effort.

By looking at the broad categories 'dependent' (FIM scores 3, 4 and 5) and 'complete dependence' (FIM scores 1 and 2), the FIM can assist in equipment needs identification (remembering that this tool is not designed for person-specific risk assessment).

NOTES

At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



ENABLE LIFECARE PTY LTD T 1300 370 370 enablelifecare.com.au