
Oxford[®] Patient Handling Belt Slings




Patient Handling Belt

The Patient Handling Belt offers carers a means of providing standing assistance to patients. It can be used independently or in combination with a compatible manual standing/transfer device, e.g. the Oxford Switch. It may also be used to help with walking rehabilitation/gait training activity following careful assessment.

- ✓ Side handles for standing assistance and support
- ✓ Adjustable belt enables use with the Oxford Switch
- ✓ Non-slip inner layer ensures secure fit
- ✓ Elasticated top band for increased comfort
- ✓ Colour coded binding indicates belt size
- ✓ Can be customised to fulfil patient specific needs



Size	Oxford	Enable	SWL
XS	SPECIAL	MTO	180 kg
S	SL3151	TBA	180 kg
M	SL3152	TBA	180 kg
L	SL3153	TBA	180 kg
XL	SPECIAL	TBA	180 kg



Chair to Chair



Bed to Chair/Chair to Bed
(From/to seated position)



Floor to Bed/Chair



Chair to Commode/Toilet



Transfer to Bath



Rehab/Assisted Standing





At Enable we source quality, affordable and innovative products globally and have them available for fast delivery. But we do more than just provide products.

We give the knowledge and support to help our customers grow their business and deliver better life care.



ENABLE LIFECARE PTY LTD

T 1300 370 370

enablelifecare.com.au