

OSKA™ Flex Series1-H

Modular foam mattress



| | |
|----------------------|-------------------|
| Prevention Treatment | High Risk |
| Comfort | N/A |
| Type | ★★★★★ Replacement |

FEATURES

- High risk rating
- Modular foam
- 4 way stretch covers
- Ideal for use on profiling beds
- Turning instructions

SPECIFICATION

| MAX USER WEIGHT | SIZE |
|-----------------|-------------|
| 213kg/33.5st | 88x198x15cm |
| Order Code | MA1010 |

'High Risk' modular foam mattress.

Series1-H enables organisations to provide an excellent 'High Risk' rated mattress for all patients. Independent studies have shown that doing this gives long term cost savings to the organisation.

The modular foam cubes react independently to conform to the patient's body, redistributing the pressure. The mattress can be turned periodically to increase its life span.

The mattress is designed to conform well on a profiling nursing bed but is also available in divan bed size to enable residential homes to provide pressure relief for all residents. This proactive approach demonstrates going above and beyond in the fight against pressure ulcers - Prevention is better than cure!

OSKA™ Flex Series1-H

Modular foam mattress

TURNING INSTRUCTIONS

January - March

Number at foot of bed: 1
To turn: Spin mattress. Do not flip.



April - June

Number at foot of bed: 2
To turn: Flip mattress and spin.



July - September

Number at foot of bed: 3
To turn: Spin mattress.



October - December

Number at foot of bed: 4
To return to original position for January:
Flip mattress and spin.

