

Mattress Selector




Involuntary movement during sleep
[approx. every 15 mins]

PATIENT MOBILITY
TISSUE DAMAGE
PRODUCT OPTIONS

Able to turn and move independently or with minimal assistance in bed.

No damage



Prevention



Fika SWL 159kg

Restricted ability
Some assisted repositioning required.

Grade 1-2 pressure damage or Grade 4 heel pressure damage




Level 1-2 Pressure Wound




Aria+ SWL 230kg

Does not move during sleep & / or spends all in bed. Minimal ability to move in bed, Assisted repositioning required on most occasions..


Grade 4 pressure damage or grade 3 pressure damage in a patient with a high/low body mass index



Level 3-4 Pressure Wound



Aria+ with Pump SWL 250kg



Promatt SWL 227kg



Alto SWL 175kg